



Column: **MEGA-5 FAST - 0.10mm, 0.10µm, 10m**  
Catalog Code: F-5-010-010-10

条件:

进样口: 分流 280°C, 1µL进样量

检测器: TOF-MS 系统

柱温箱程序: 50-300°C (2min), 35°C/min

载气: 氦气, 0.6mL/min

1	4.360 min	Cotinine	可替宁
2	4.783 min	Caffeine	咖啡因
3	6.153 min	SKF	
4	6.480 min	Codeine (TMS)	可待因
5	6.607 min	Morphine (TMS)	吗啡
6	6.753 min	6-MAM (TMS)	单乙酰吗啡
7	6.887 min	Nalorphine (TMS)	烯丙吗啡

Analysis carried out with a DANI Master GC coupled with DANI Master TOF MS detector.