

# 溶剂 on MEGA-SOLVE Columns

Column: **MEGA-SOLVE 1 - 0.32mm, 0.25 $\mu$ m, 50m**

Catalog Code: C-SOLVE1-032-025-50

**MEGA-SOLVE Columns are especially developed for analysis of very complex mixtures of solvents.**

条件:

进样口: 分流 200 $^{\circ}$ C, 1:100分流比, 0.5 $\mu$ L进样量

检测器: FID 200 $^{\circ}$ C

柱温箱程序: 35-90 $^{\circ}$ C, 1 $^{\circ}$ C/min

载气: 氦气, 2.5mL/min

**Peaks:**

- |                         |                       |                               |
|-------------------------|-----------------------|-------------------------------|
| 1. Ethyl Ether          | 8. Trichloro Ethylene | 17. p-Xylene                  |
| 2. Acetone              | 9. Ethanol            | 18. Mesityl Oxide             |
| 3. Carbon Tetrachloride | 10. Acetonitrile      | 19. o-Xylene                  |
| 4. Ethyl Acetate        | 11. Toluene           | 20. n-Butanol                 |
| 5. Methylene Chloride   | 12. 1,4-dioxane       | 21. Pyridine                  |
| 6. Benzene              | 13. n-Butyl Acetate   | 22. 2-Methoxy-Ethanol         |
| 7. Methanol             | 14. Propanol          | 23. 2-Methoxy-Ethanol Acetate |
|                         | 15. Ethyl Benzene     | 24. 2-Ethoxy-Ethanol          |
|                         | 16. m-Xylene          | 25. 2-Ethoxy-Ethanol Acetate  |
|                         |                       | 26. N,N-dimethylformamide     |
|                         |                       | 27. Diacetone Alcohol         |
|                         |                       | 28. N,N-dimethyl Aniline      |
|                         |                       | 29. N,N-dimethyl Acetamide    |
|                         |                       | 30. N,N-diethyl Aniline       |

