



溶剂

on MEGA-VOC 1 column (USP G43 phase)

Column: **MEGA-VOC 1 - 0.25mm, 1.50 μ m, 30m**

Catalog Code: C-VOC1-025-150-30

条件:

进样口: 分流 250 $^{\circ}$ C, 150mL/min分流比, 0.5 μ L进样量

检测器: FID 290 $^{\circ}$ C

柱温箱程序: 40-260 $^{\circ}$ C, 1 $^{\circ}$ C/min

载气: 氦气, 0.7mL/min

Peaks:

1. Methanol
2. Ethanol
3. Acetone
4. Isopropyl Alcohol
5. Dichloromethane
6. n-Hexane
7. MEK (Methyl Ethyl Ketone)
8. Ethyl Acetate
9. Cyclohexane
10. 1,4-dioxane
11. Toluene
12. Cyclopentanone
13. N-methylpyrrolidone